

BEDS

Purpose

This document provides the weight rating when actuating the Schwintek Bed Tilt System.

System Rating

The Schwintek Bed Tilt system is designed and rated to be actuated with no more than **800 lbs (363 Kgs)** of total weight. Using the system with weights higher than recommended will void the warranty and will cause damage to the system.

Bed Tilt Lowered Position



Bed Tilt Raised Position



As a supplier of components to the RV industry, safety, education and customer satisfaction are our primary concerns. Should you have any questions, please do not hesitate to contact us at (574) 537-8900 or by email at customerservice@lci1.com. Self-help tips, technical documents, product videos and a training class schedule are available at lippert.com or by downloading the MyLCI app.



This manual is provided courtesy of **RV Tech Cheat Sheet**



Professional resources for the whole RV industry — free.

We put manuals, wiring diagrams and hard-won tech knowledge in every technician's hands, so repairs get done right the first time, techs work faster, and RV owners everywhere get better service.



rvtechcheatsheet.com/manuals

Find thousands more RV manuals, wiring diagrams and parts lists

RV techs: share fixes, ask questions and learn from each other.

Join our free Facebook group →



Follow proper safety practices. If in doubt, consult a professional RV technician.